



### CF-400 Precision Co-Efficient of Friction Tester

The cf-400 precision co-efficient of friction tester measures the static and kinetic friction properties of plastic films, foils, laminates, papers and boards. The equipment performs tests to recognised international test standards including Astm d1894, iso 8295, bs 2782 method 824a, and tappi t549.

This equipment is essential for measuring the slip properties of packaging materials to ensure smooth running on production packaging machines or to measure the effect that a coating or print has on base material.

The CF-400 features the latest in design and technology for simple set up, testing, measurement and recording using a touchscreen display and data output via USB.

The constant smooth lead screw driven cross arm, standard base material with controlled surface roughness, and precision load cell ensures reliable and repeatable measurement. Test results are displayed on the touchscreen for Static and Dynamic Friction, Active Force, Peak Force and Minimum Force in grams to 3 decimal places, during and after the test a test graph can be viewed.

Test data is automatically saved with date/time stamp for download to .csv for further analysis. A standard Excel template can be used to process data for test reports.

### Other benefits include:

Quick action Spring Clamp for base sample material (optional vacuum channel clamp), easy fit sled with magnetic clamps for sample material, precision load cell with link connectors, standard base material with controlled surface roughness, optional replacement test beds in various materials, optional accessories for peel, tear and tensile tests.



### UK & INTERNATIONAL

RDM Test Equipment Ltd.  
39 Gold Nurseries Business Park  
Jenkins Drive  
Elsternham  
Hertfordshire  
CM22 6JX  
Tel: +44 (0)1279 817171  
Email: sales@rdmtest.com  
Website: rdmtest.com

### AMERICAS

RDM Test Equipment LLC.  
441 Old Hwy 8 NW  
Suite 203  
New Brighton  
Minneapolis  
MN55112  
USA  
Tel: +1 651 766 2565  
Mobile: +1 651 485 2372  
Email: sales@rdmtest.com  
Website: rdmtest.com

## Specification

**Test Bed:** Natural anodised cast aluminium, 120mm x 420mm with maximum 150mm travel. Spring clamp fitted as standard, interchangeable with calibration pulley and optional peel/tear/tensile grips.

**Test Area:** 374 x 100mm

**Test Sled:** Anodised aluminium with foam contact pad. 200grams, 63.5 mm x 63.5mm (2.5" x 2.5"). Solid and cord links. Other sleds available.

**Test Arm:** Fitted with precision load cell, 2kg capacity. Link pin fitting interchangeable with optional peel. tear, tensile grips.

**Speed Control:**  
10 – 1000mm/min + / - 1% FS

**Motor Drive:** DC synchronous motor driving ball screw and test arm.

**Force:** 0-1000 grammes + / - 0.25% FS. Active display on touchscreen.

**Distance:** 0 to 150mm + / - 1% FS. Active display on touchscreen.

**Test Results:** Static & Dynamic Friction (COF) calculated value 0.000 to 1.000 + / - 0.25% FS. Peak Force, Minimum Force.

**Touch Panel Screen:**  
7" Touchscreen display Kinco HMI  
16.77m (24 bits) True Colour 32-bit 792 Mhz Arm Risc CPU.

## Optional Vacuum Test Bed:

If fitted, requires compressed air input pressure of 80 – 100 PSI, supplied to vacuum generator with on/off push button. Vacuum drawn on test bed to hold base sample.

**Test Bed Temperature:** Ambient

**Dimensions:**  
440 x 265 x 275mm approx. 25kg

**Data Output:** Data formatted to .csv output to USB

**Power:** 240 VAC single phase 50/60 Hz (110V AC available on request) 0.75 KW max

**Environment:** 5-35°C ambient operating temperature, RH 75% max (non- condensing)

## Standard Accessories

Test Sled, 200g, 63.5mm x 63.5mm.  
Sled / Bed Cutting Templates.  
Check weight & pulley.  
Sled links.

## Options

Ski sled 100g for measurement of stainless steel on test material.  
Bed replacements to give test comparison with different metals.  
Tensile grips for T peel, tear or tensile test.  
Peel attachment for 90° and 180° label / sticky tape peel test.



**UK & INTERNATIONAL**  
RDM Test Equipment Ltd.  
39 Gold Nurseries Business Park  
Jenkins Drive  
Elsenham  
Hertfordshire  
CM22 6JX  
Tel: +44 (0)1279 817171  
Email: sales@rdmtest.com  
Website: rdmttest.com

**AMERICAS**  
RDM Test Equipment LLC.  
441 Old Hwy 8 NW  
Suite 203  
New Brighton  
Minneapolis  
MN55112  
USA  
Tel: +1 651 766 2565  
Mobile: +1 651 485 2372  
Email: sales@rdmttest.com  
Website: rdmttest.com